


INSTRUMENTS DOLORIMETER



Sto	Contents ock kit #
	1.52cm² flat circular probe
0	2 pound (.90 kg) gauge #12-1440 5 pound gauge (2.27 kg) #12-1391 25 pound gauge (11.34 kg) #12-1442 60 pound gauge (27.22 kg) #12-1143

Baseline® Dolorimeter

Evaluate p ain threshold. S tandardized (1.52cm ²) flat circular probe is pushed against subject's skin until pain threshold is reached. Gauge reads in pounds and kilograms.